

Impact on gas-permeable contact lens parameters after storage in a non-neutralized hydrogen peroxide case during 1 to 30 days.

Dr Langis Michaud O.D. M.Sc. FAAO (Dipl) FSLs FBCLA FEAOO
Professor



INTRODUCTION

Hydrogene peroxide is more widely use to clean and store rigid gas permeable lenses

Scleral lenses are prescribed to treat eye disease. Patients may benefit from a hydrogen peroxide-based preservative-free care regimen. However, the size of the lenses does not fit into the basket provided by the manufacturer.

Logically, because scleral lenses does not absorb solution, a regular contact lens case may be used, provided that a thoroughful rinsing step is done before lens wear.

There were few studies on the impact on contact lens parameters of long-term in a non-neutralized peroxide solution.




PURPOSE

To evaluate the impact of a non-neutralized hydrogen peroxide solution on scleral lens parameters

METHODS

Clean diagnostic scleral lenses (Hexafocon B) were stored for 1,3,7, and 30 days in a regular contact lens case, with a hydrogen peroxide solution (Clear Care, Alcon US)



Lens parameters were assessed after lens soaking (base curve, power, diameter)

Lens wetting angle was measured – sessile drop technique

Hydrogene Peroxide concentration was measured.


RESULTS

No change in lens parameters during the experimentation


Wetting angle improved by Day 7

$F(4,72) = < 0,001; p < 0,05; \omega^2 = 0,448$

Hydrogene Peroxide concentration remained stable over time



Wetting angle



H2O2 concentration

CONCLUSION

Hydrogene Peroxide can be used in regular non -neutralizing case with no impact on lens parameters

May represent an interesting option for patients wearing larger scleral lenses . May represent a very good alternative for lens storage in office